

Package ‘mSigHdp’

October 30, 2020

Title Mutational signature extraction using hdp (Hierarchical Dirichlet Process)

Version 1.1.2001

Description Adapts HDP to mutational signature extraction. Calls hdp for mutational signature analysis. mSigHdp incorporates burnin and Gibbs sampling processes, followed up by extracting signatures from raw clusters collected by Gibbs sampling, and improved diagnostic plotting with visualization on signature and exposures.

License GPL-3

Encoding UTF-8

LazyData true

Language en-US

BuildManual no

biocViews

Roxygen list(markdown = TRUE)

Depends R (>= 4.0)

RoxygenNote 7.1.1

Remotes github::steverozen/hdp@*release,
github::steverozen/ICAMSxtra@*release

Imports hdp (>= 0.3.0),
ICAMS (>= 2.2.4),
reshape2

Suggests ICAMSxtra (>= 0.0.2),
testthat,
utils

R topics documented:

AnalyzeAndPlotretval	2
ChainBurnin	3
ChainsDiagnosticPlot	4
CleanChlist	4
CombineChainsAndExtractSigs	5
CombinePosteriorChains	6
ComponentDiagnosticPlotting	8

ExtendBurnin	9
GenerateAverageCluster	9
Generateppindex	10
GeneratePriorppindex	10
MultipleSetupAndPosterior	11
OldRunHdpParallel	13
PlotSamplesHighSigExp	17
PrepInit	18
PriorSetupAndActivate	19
RunHdpxParallel	20
SetupAndActivate	24
SetupAndPosterior	25
test.ground.truth.sig	27
test.spectra	28

Index 29

AnalyzeAndPlotretval

*Evaluate and plot retval from CombinePosteriorChains or
CombineChainsAndExtractSigs This function now works for
both NR's pipeline and Mo's pipeline*

Description

Evaluate and plot retval from `CombinePosteriorChains` or `CombineChainsAndExtractSigs`
This function now works for both NR's pipeline and Mo's pipeline

Usage

```
AnalyzeAndPlotretval(
  retval,
  input.catalog,
  out.dir = NULL,
  ground.truth.sig = NULL,
  ground.truth.exp = NULL,
  verbose = TRUE,
  overwrite = TRUE,
  diagnostic.plot = TRUE
)
```

Arguments

<code>retval</code>	the output from function <code>CombinePosteriorChains</code> or <code>CombineChainsAndExtractSigs</code>
<code>input.catalog</code>	input catalog matrix or path to file with input catalog
<code>out.dir</code>	Directory that will be created for the output; if <code>overwrite</code> is <code>FALSE</code> then abort if <code>out.dir</code> already exists.
<code>ground.truth.sig</code>	Optional. Either a string with the path to file with ground truth signatures or and ICAMS catalog with the ground truth signatures. These are the signatures used to construct the ground truth spectra.

<code>ground.truth.exp</code>	Optional. Ground truth exposure matrix or path to file with ground truth exposures. If NULL skip checks that need this information.
<code>verbose</code>	If TRUE then message progress information.
<code>overwrite</code>	If TRUE overwrite <code>out.dir</code> if it exists, otherwise raise an error.
<code>diagnostic.plot</code>	If TRUE plot diagnostic plot. This is optional because there are cases having error

ChainBurnin	<i>Prepare an hdpState-class object and run the Gibbs sampling burnin.</i>
-------------	--

Description

Prepare an [hdpState-class](#) object and run the Gibbs sampling burnin.

Usage

```
ChainBurnin(
  hdp.state,
  seedNumber = 1,
  burnin = 5000,
  cpiter = 3,
  burnin.verbosity = 0,
  burnin.multiplier = 2,
  burnin.checkpoint = TRUE
)
```

Arguments

<code>hdp.state</code>	An hdpState-class object or a list representation of an hdpState-class object.
<code>seedNumber</code>	An integer that is used to generate separate random seeds for the call to dp_activate , and before the call of hdp_burnin .
<code>burnin</code>	Pass to hdp_burnin burnin. The number of burn-in iterations
<code>cpiter</code>	Pass to hdp_burnin cpiter. The number of iterations of concentration parameter sampling to perform after each iteration.
<code>burnin.verbosity</code>	Pass to hdp_burnin verbosity. Verbosity of debugging statements.
<code>burnin.multiplier</code>	A checkpoint setting. <code>burnin.multiplier</code> rounds of burnin iterations will be run. After each round, a burn-in chain will be save for checkpoint. A total number of 10,000 iterations is recommended for most analysis. However, number of iterations can be adjusted based on the size of dataset. According to our experience, 20,000 iterations are needed when analyzing all PCAWG7 genomes (2,780 samples). The burnin can be continued from a checkpoint file with ExtendBurnin .
<code>burnin.checkpoint</code>	Default is False. If True, a checkpoint for burnin will be created.

Value

A list with 2 elements:

`hdplist` A list representation of an `hdpState-class` object.

likelihood A numeric vector with the likelihood at each iteration.

`ChainsDiagnosticPlot`

Diagnostic plot for a `hdpSampleMulti` object

Description

Diagnostic plot for a `hdpSampleMulti` object

Usage

```
ChainsDiagnosticPlot(retval, input.catalog, out.dir, verbose)
```

Arguments

<code>retval</code>	<p>output from <code>CombinePosteriorChains</code>. A list with the following elements:</p> <p>signature The extracted signature profiles as a matrix; rows are mutation types, columns are samples (e.g. tumors).</p> <p>exposure The inferred exposures as a matrix of mutation counts; rows are signatures, columns are samples (e.g. tumors).</p> <p>multi.chains A <code>hdpSampleMulti-class</code> object. This object has the method <code>chains</code> which returns a list of <code>hdpSampleChain-class</code> objects. Each of these sample chains objects has a method <code>final_hdpState</code> (actually the methods seems to be just <code>hdp</code>) that returns the <code>hdpState</code> from which it was generated.</p>
<code>input.catalog</code>	ground truth catalog
<code>out.dir</code>	Directory that will be created for the output; if <code>overwrite</code> is <code>FALSE</code> then abort if <code>out.dir</code> already exists.
<code>verbose</code>	If <code>TRUE</code> then message progress information.

`CleanChlist`

If the job of Gibbs sampling from `MultipleSetupAndPosterior` has an error caught by R, the corresponding element of `chlist` has class `try-error`. If the job is stopped with, e.g. a segfault, the `chlist` element is `NULL`.

Description

If the job of Gibbs sampling from `MultipleSetupAndPosterior` has an error caught by R, the corresponding element of `chlist` has class `try-error`. If the job is stopped with, e.g. a segfault, the `chlist` element is `NULL`.

Usage

```
CleanChlist(chlist, verbose = FALSE)
```

Arguments

`chlist` A list of `hdpSampleChain-class` objects. Return from `MultipleSetupAndPosterior`.
`verbose` If TRUE then message progress information.

Value

Invisibly, the clean, non-error `chlist` This is a list of `hdpSampleChain-class` objects.

```
CombineChainsAndExtractSigs
```

Extract components and exposures from multiple posterior sample chains This function returns signatures with high confidence (found in more than 90% #' posterior samples)

Description

Extract components and exposures from multiple posterior sample chains This function returns signatures with high confidence (found in more than 90% #' posterior samples)

Usage

```
CombineChainsAndExtractSigs(
  clean.chlist,
  input.catalog,
  multi.types,
  verbose = TRUE,
  cos.merge = 0.9,
  confident.prop = 0.9,
  noise.prop = 0.1,
  hc.cutoff = 0.1
)
```

Arguments

`clean.chlist` A list of `hdpSampleChain-class` objects. Each element is the result of one posterior sample chain.
`input.catalog` Input spectra catalog as a matrix or in `ICAMS` format.
`multi.types` A logical scalar or a character vector. If FALSE, The HDP analysis will regard all input spectra as one tumor type. HDP structure as one parent node for all tumors
 If TRUE, the HDP analysis will infer tumor types based on the string before "::" in their names. e.g. tumor type for "SA.Syn.Ovary-AdenoCA::S.500" would be "SA.Syn.Ovary-AdenoCA" HDP structure as a grandparent node for whole data and one parent node for each tumor type

	<p>If <code>multi.types</code> is a character vector, then it should be of the same length as the number of columns in <code>input.catalog</code>, and each value is the name of the tumor type of the corresponding column in <code>input.catalog</code>. e.g. <code>c("SA.Syn.Ovary-AdenoCA", "SA.Syn.Kidney-RCC")</code>.</p>
<code>verbose</code>	If TRUE then message progress information.
<code>cos.merge</code>	The cosine similarity threshold for merging raw clusters from the posterior sampling chains into "components" i.e. signatures; passed to extract_components_from_clusters
<code>confident.prop</code>	Pass to interpret_components . clusters with at least <code>confident.prop</code> of posterior samples are high confident signatures
<code>noise.prop</code>	Pass to interpret_components . Clusters with less than <code>noise.prop</code> of posterior samples are noise signatures
<code>hc.cutoff</code>	Pass to extract_components_from_clusters . The cutoff of height of hierarchical clustering dendrogram

Value

Invisibly, a list with the following elements:

- signature** The extracted signature profiles as a matrix; rows are mutation types, columns are signatures including high confident signatures - 'hdp' signatures and moderate confident signatures - 'potential hdp' signatures.
- signature.post.samp.number** A dataframe with two columns. The first column corresponds to each signature in `signature` and the second columns contains the number of posterior samples that found the raw clusters contributing to the signature.
- signature.cdc** A [comp_dp_counts](#) like dataframe. Each column corresponds to the sum of all [comp_dp_counts](#) matrices of the raw clusters contributing to each signature in `signature`
- exposureProbs** The inferred exposures as a matrix of mutation probabilities; rows are signatures, columns are samples (e.g. tumors).
- noise.signature** The extracted signature profiles as a matrix; rows are mutation types, columns are signatures with less than `noise.prop` of posterior samples.
- noise.post.samp.number** A data frame with two columns. The first column corresponds to each signature in `noise.signature` and the second columns contains the number of posterior samples that found the raw clusters contributing to the signature.
- noise.cdc** A [comp_dp_counts](#) like data frame. Each column corresponds to the sum of all [comp_dp_counts](#) matrices of the raw clusters contributing to each signature in `code-noise.signature`
- extracted.retval** A list object returned from code [interpret_components](#).

CombinePosteriorChains

Extract components and exposures from multiple posterior sample chains

Description

Extract components and exposures from multiple posterior sample chains

Usage

```
CombinePosteriorChains(
  clean.chlist,
  input.catalog,
  multi.types,
  verbose = TRUE,
  cos.merge = 0.9,
  categ.CI = 0.95,
  exposure.CI = 0.95,
  min.sample = 1,
  diagnostic.folder = NULL
)
```

Arguments

- `clean.chlist` A list of [hdpSampleChain-class](#) objects. Each element is the result of one posterior sample chain.
- `input.catalog` Input spectra catalog as a matrix or in [ICAMS](#) format.
- `multi.types` A logical scalar or a character vector. If `FALSE`, The HDP analysis will regard all input spectra as one tumor type.
If `TRUE`, the HDP analysis will infer tumor types based on the string before ":" in their names. e.g. tumor type for "SA.Syn.Ovary-AdenoCA::S.500" would be "SA.Syn.Ovary-AdenoCA"
If `multi.types` is a character vector, then it should be of the same length as the number of columns in `input.catalog`, and each value is the name of the tumor type of the corresponding column in `input.catalog`.
e.g. `c("SA.Syn.Ovary-AdenoCA", "SA.Syn.Kidney-RCC")`.
- `verbose` If `TRUE` then message progress information.
- `cos.merge` The cosine similarity threshold for merging raw clusters from the posterior sampling chains into "components" i.e. signatures; passed to [hdp_extract_components](#).
- `categ.CI` A number the range `[0, 1]`. The level of the confidence interval used in step 4 of [hdp_merge_and_extract_components](#). This governs when "averaged raw cluster" get assigned to component 0, i.e. if the the confidence interval overlaps 0. Lower values make it less likely that an averaged raw cluster will be assigned to component 0. The CI in question is for the number of mutations in a given mutation class (e.g. `ACA > AAA`, internally called a "category"). If, for every mutation class, this CI overlaps 0, then the averaged raw cluster goes to component 0.
- `exposure.CI` A number in the range `[0, 1]`. The level of the confidence interval used in step 5 of [hdp_merge_and_extract_components](#). The CI in question here for the total number of mutations assigned to an averaged raw cluster.
- `min.sample` A "component" (i.e. signature) must have at least this many samples; passed to [hdp_merge_and_extract_components](#).
- `diagnostic.folder` If provided, diagnostic plots for `hdp.0` components are provided

Value

Invisibly, a list with the following elements:

signature The extracted signature profiles as a matrix; rows are mutation types, columns are samples (e.g. tumors).

exposure The inferred exposures as a matrix of mutation counts; rows are signatures, columns are samples (e.g. tumors).

multi.chains A `hdpSampleMulti-class` object. This object has the method `chains` which returns a list of `hdpSampleChain-class` objects. Each of these sample chains objects has a method `final_hdpState` (actually the methods seems to be just `hdp`) that returns the `hdpState` from which it was generated.

sum_raw_clusters_after_cos_merge A matrix containing aggregated spectra of raw clusters after cosine similarity merge step in `hdp_merge_and_extract_components`.

sum_raw_clusters_after_nonzero_categ A matrix containing aggregated spectra of raw clusters after non-zero category selecting step in `hdp_merge_and_extract_components`.

clust_hdp0_ccc4 A matrix containing aggregated spectra of raw clusters moving to `hdp.0` after non-zero category selection step in `hdp_merge_and_extract_components`.

clust_hdp0_ccc5 A matrix containing aggregated spectra of raw clusters moving to `hdp.0` after non-zero observation selection step in `hdp_merge_and_extract_components`.

ComponentDiagnosticPlotting

Diagnostic plot for a `hdpSampleMulti` object. This function is compatible with the return object from Liu's `extract_components_from_clusters`

Description

Diagnostic plot for a `hdpSampleMulti` object. This function is compatible with the return object from Liu's `extract_components_from_clusters`

Usage

```
ComponentDiagnosticPlotting(retval, input.catalog, out.dir, verbose)
```

Arguments

<code>retval</code>	Return from <code>CombineChainsAndExtractSigs</code>
<code>input.catalog</code>	Input spectra catalog as a matrix or in <code>ICAMS</code> format.
<code>out.dir</code>	Directory that will be created for the output; if <code>overwrite</code> is <code>FALSE</code> then abort if <code>out.dir</code> already exists.
<code>verbose</code>	If <code>TRUE</code> then message progress information.

ExtendBurnin	<i>Extend Burn in iteration for a list representation of an <code>hdpState-class</code> object. This list is an output from <code>hdp_burnin</code> or <code>ActivateandBurnin</code>.</i>
--------------	--

Description

Extend Burn in iteration for a list representation of an `hdpState-class` object. This list is an output from `hdp_burnin` or `ActivateandBurnin`.

Usage

```
ExtendBurnin(hdplist, seedNumber = 1, burnin = 4000, cpiter = 3, verbosity = 0)
```

Arguments

<code>hdplist</code>	A list representation of an <code>hdpState-class</code> object
<code>seedNumber</code>	A random seed for setting the environment of <code>hdp_burnin</code> .
<code>burnin</code>	Pass to <code>hdp_posterior</code> burnin.
<code>cpiter</code>	Pass to <code>hdp_posterior</code> cpiter.
<code>verbosity</code>	Pass to <code>hdp_posterior</code> verbosity.

Value

A list with hdp object after burn-in iteration and likelihood of iteration

GenerateAverageCluster	<i>Generate average pattern of clusters of each posterior chain from combined list of multiple posterior sample chains</i>
------------------------	--

Description

Generate average pattern of clusters of each posterior chain from combined list of multiple posterior sample chains

Usage

```
GenerateAverageCluster(clean.chlist)
```

Arguments

<code>clean.chlist</code>	A list of multiple (or one) posterior sample chains.
---------------------------	--

Value

A list of matrices containing the average pattern of clusters within each posterior chain and a list of matrices containing the sum of each cluster in each posterior chain

Generateppindex	<i>Generate index for a HDP structure and num.tumor.types for other functions</i>
-----------------	---

Description

Generate index for a HDP structure and num.tumor.types for other functions

Usage

```
Generateppindex(multi.types, input.catalog)
```

Arguments

`multi.types` A logical scalar or a character vector. If FALSE, The HDP analysis will regard all input spectra as one tumor type. HDP structure as one parent node for all tumors

If TRUE, the HDP analysis will infer tumor types based on the string before ":" in their names. e.g. tumor type for "Ovary-AdenoCA::S.500" would be "Ovary-AdenoCA" HDP structure as a grandparent node for whole data and one parent node for each tumor type

If `multi.types` is a character vector, then it should be of the same length as the number of columns in `input.catalog`, and each value is the name of the tumor type of the corresponding column in `input.catalog`.

If not FALSE, HDP will give a parent node for each tumor type and a grandparent node for the whole dataset.

e.g. `c("Ovary-AdenoCA", "Kidney-RCC")`.

`input.catalog`

Input spectra catalog as a matrix or in [ICAMS](#) format.

GeneratePriorppindex	<i>Generate index for a HDP structure and num.tumor.types for other functions for hdp_prior_init</i>
----------------------	--

Description

Generate index for a HDP structure and num.tumor.types for other functions for hdp_prior_init

Usage

```
GeneratePriorppindex(multi.types, input.catalog, nps)
```

Arguments

<code>multi.types</code>	<p>A logical scalar or a character vector. If <code>FALSE</code>, The HDP analysis will regard all input spectra as one tumor type.</p> <p>If <code>TRUE</code>, the HDP analysis will infer tumor types based on the string before ":" in their names. e.g. tumor type for "SA.Syn.Ovary-AdenoCA::S.500" would be "SA.Syn.Ovary-AdenoCA"</p> <p>If <code>multi.types</code> is a character vector, then it should be of the same length as the number of columns in <code>input.catalog</code>, and each value is the name of the tumor type of the corresponding column in <code>input.catalog</code>. e.g. <code>c("SA.Syn.Ovary-AdenoCA", "SA.Syn.Kidney-RCC")</code>.</p>
<code>input.catalog</code>	Input spectra catalog as a matrix or in ICAMS format.
<code>nps</code>	Number of prior signatures

MultipleSetupAndPosterior

Activate hierarchical Dirichlet processes and run posterior sampling in parallel.

Description

Activate hierarchical Dirichlet processes and run posterior sampling in parallel.

Usage

```
MultipleSetupAndPosterior(
  input.catalog,
  seedNumber = 1,
  K.guess,
  multi.types = FALSE,
  verbose = TRUE,
  burnin = 5000,
  burnin.multiplier = 2,
  burnin.checkpoint = TRUE,
  post.n = 200,
  post.space = 100,
  post.cpiter = 3,
  post.verbosity = 0,
  CPU.cores = 20,
  num.child.process = 20,
  gamma.alpha = 1,
  gamma.beta = 20,
  gamma0.alpha = gamma.alpha,
  gamma0.beta = gamma.beta,
  checkpoint.chlist = TRUE,
  checkpoint.l.chain = TRUE,
  prior.sigs = NULL,
  prior.pseudoc = NULL
)
```

Arguments

<code>input.catalog</code>	Input spectra catalog as a matrix or in ICAMS format.
<code>seedNumber</code>	A random seeds passed to dp_activate .
<code>K.guess</code>	Suggested initial value of the number of clusters. Usually, the number of clusters is two times of the number of extracted signatures. Passed to dp_activate as <code>initcc</code> .
<code>multi.types</code>	<p>A logical scalar or a character vector. If <code>FALSE</code>, The HDP analysis will regard all input spectra as one tumor type. HDP structure as one parent node for all tumors</p> <p>If <code>TRUE</code>, the HDP analysis will infer tumor types based on the string before "::" in their names. e.g. tumor type for "Ovary-AdenoCA::S.500" would be "Ovary-AdenoCA" HDP structure as a grandparent node for whole data and one parent node for each tumor type</p> <p>If <code>multi.types</code> is a character vector, then it should be of the same length as the number of columns in <code>input.catalog</code>, and each value is the name of the tumor type of the corresponding column in <code>input.catalog</code>.</p> <p>If not <code>FALSE</code>, HDP will give a parent node for each tumor type and a grandparent node for the whole dataset.</p> <p>e.g. <code>c("Ovary-AdenoCA", "Kidney-RCC")</code>.</p>
<code>verbose</code>	If <code>TRUE</code> then message progress information.
<code>burnin</code>	Pass to hdp_burnin <code>burnin</code> . The number of burn-in iterations
<code>burnin.multiplier</code>	A checkpoint setting. <code>burnin.multiplier</code> rounds of <code>burnin</code> iterations will be run. After each round, a burn-in chain will be save for checkpoint. A total number of 10,000 iterations is recommended for most analysis. However, number of iterations can be adjusted based on the size of dataset. According to our experience, 20,000 iterations are needed when analyzing all PCAWG7 genomes (2,780 samples). The burnin can be continued from a checkpoint file with ExtendBurnin .
<code>burnin.checkpoint</code>	Default is <code>False</code> . If <code>True</code> , a checkpoint for burnin will be created.
<code>post.n</code>	Pass to hdp_posterior_sample <code>n</code> . The number of posterior samples to collect.
<code>post.space</code>	Pass to hdp_posterior_sample <code>space</code> . The number of iterations between collected samples.
<code>post.cpiter</code>	Pass to hdp_posterior_sample and hdp_burnin <code>cpiter</code> . The number of iterations of concentration parameter sampling to perform after each iteration
<code>post.verbosity</code>	Pass to hdp_posterior_sample <code>verbosity</code> . Verbosity of debugging statements. No need to change unless for development purpose
<code>CPU.cores</code>	Number of CPUs to use; this should be no more than <code>num.child.process</code> .
<code>num.child.process</code>	Number of posterior sampling chains; can set to 1 for testing. We recommend 20 for real data analysis
<code>gamma.alpha</code>	Shape parameter of the gamma distribution prior for the Dirichlet process concentration parameters; in this function the gamma distributions for all Dirichlet processes, except possibly the top level process, are the same.

<code>gamma.beta</code>	<p>Inverse scale parameter (rate parameter) of the gamma distribution prior for the Dirichlet process concentration parameters; in this function the gamma distributions for all Dirichlet processes, except possibly the top level process, are the same.</p> <p>We recommend <code>gamma.alpha = 1</code> and <code>gamma.beta = 20</code> for single-base-substitution signatures extraction; <code>gamma.alpha = 1</code> and <code>gamma.beta = 50</code> for doublet-base-substitution/INDEL signature extraction</p>
<code>gamma0.alpha</code>	See figure B.1 from Nicola Robert's thesis. The shape parameter (α_0) of the gamma distribution priors for the Dirichlet process concentration parameters (γ_0) for G_0 .
<code>gamma0.beta</code>	See figure B.1 from Nicola Robert's thesis. Inverse scale parameter (rate parameter, β_0) of the gamma distribution priors for the Dirichlet process concentration parameters (γ_0) for G_0 .
<code>checkpoint.chlist</code>	If TRUE, checkpoint the (unclean) chlist to "initial.chlist.Rdata" in the current working directory. and checkpoint the clean chlist to "clean.chlist.Rdata" in the current working directory.
<code>checkpoint.1.chain</code>	If TRUE checkpoint the sample chain to current working directory, in a file called <code>sample.chain.seed_number.Rdata</code> .
<code>prior.sigs</code>	A matrix containing prior signatures.
<code>prior.pseudoc</code>	<p>A numeric list. Pseudo counts of each prior signature. Recommended is 1000. In practice, it may be advisable to put lower weights on prior signatures that you do not expect to be present in your dataset, or even exclude some priors entirely.</p>

Value

Invisibly, the clean `chlist` (output of `CleanChlist`). This is a list of `hdpSampleChain-class` objects.

`OldRunHdpParallel` *Deprecated, extract mutational signatures and optionally compare them to existing signatures and exposures.*

Description

Deprecated, This functions uses the original method of combining raw clusters into "components". Use `RunHdpxParallel` instead.

Usage

```
OldRunHdpParallel(
  input.catalog,
  seedNumber = 1,
  K.guess,
  multi.types = FALSE,
  verbose = TRUE,
  post.burnin = 5000,
```

```

post.n = 200,
post.space = 100,
post.cputer = 3,
post.verbosity = 0,
CPU.cores = 20,
num.child.process = 20,
cos.merge = 0.9,
min.sample = 1,
categ.CI = 0.95,
exposure.CI = 0.95,
ground.truth.sig = NULL,
ground.truth.exp = NULL,
overwrite = TRUE,
out.dir = NULL,
gamma.alpha = 1,
gamma.beta = 20,
gamma0.alpha = gamma.alpha,
gamma0.beta = gamma.beta,
checkpoint.chlist = TRUE,
checkpoint.l.chain = TRUE,
prior.sigs = NULL,
prior.pseudoc = NULL,
burnin.multiplier = 2,
burnin.checkpoint = TRUE
)

```

Arguments

<code>input.catalog</code>	Input spectra catalog as a matrix or in ICAMS format.
<code>seedNumber</code>	A random seeds passed to dp_activate .
<code>K.guess</code>	Suggested initial value of the number of clusters. Usually, the number of clusters is two times of the number of extracted signatures. Passed to dp_activate as <code>initcc</code> .
<code>multi.types</code>	<p>A logical scalar or a character vector. If <code>FALSE</code>, The HDP analysis will regard all input spectra as one tumor type. HDP structure as one parent node for all tumors</p> <p>If <code>TRUE</code>, the HDP analysis will infer tumor types based on the string before ":" in their names. e.g. tumor type for "Ovary-AdenoCA::S.500" would be "Ovary-AdenoCA" HDP structure as a grandparent node for whole data and one parent node for each tumor type</p> <p>If <code>multi.types</code> is a character vector, then it should be of the same length as the number of columns in <code>input.catalog</code>, and each value is the name of the tumor type of the corresponding column in <code>input.catalog</code>.</p> <p>If not <code>FALSE</code>, HDP will give a parent node for each tumor type and a grandparent node for the whole dataset.</p> <p>e.g. <code>c("Ovary-AdenoCA", "Kidney-RCC")</code>.</p>
<code>verbose</code>	If <code>TRUE</code> then message progress information.
<code>post.n</code>	Pass to hdp_posterior_sample <code>n</code> . The number of posterior samples to collect.

<code>post.space</code>	Pass to <code>hdp_posterior_sample</code> space. The number of iterations between collected samples.
<code>post.cptiter</code>	Pass to <code>hdp_posterior_sample</code> and <code>hdp_burnin</code> <code>cptiter</code> . The number of iterations of concentration parameter sampling to perform after each iteration
<code>post.verbosity</code>	Pass to <code>hdp_posterior_sample</code> <code>verbosity</code> . Verbosity of debugging statements. No need to change unless for development purpose
<code>CPU.cores</code>	Number of CPUs to use; this should be no more than <code>num.child.process</code> .
<code>num.child.process</code>	Number of posterior sampling chains; can set to 1 for testing. We recommend 20 for real data analysis
<code>cos.merge</code>	The cosine similarity threshold for merging raw clusters from the posterior sampling chains into "components" i.e. signatures; passed to <code>hdp_extract_components</code> .
<code>min.sample</code>	A "component" (i.e. signature) must have at least this many samples; passed to <code>hdp_merge_and_extract_components</code> .
<code>categ.CI</code>	A number the range [0, 1]. The level of the confidence interval used in step 4 of <code>hdp_merge_and_extract_components</code> . This governs when "averaged raw cluster" get assigned to component 0, i.e. if the the confidence interval overlaps 0. Lower values make it less likely that an averaged raw cluster will be assigned to component 0. The CI in question is for the number of mutations in a given mutation class (e.g. ACA > AAA, internally called a "category"). If, for every mutation class, this CI overlaps 0, then the averaged raw cluster goes to component 0.
<code>exposure.CI</code>	A number in the range [0, 1]. The level of the confidence interval used in step 5 of <code>hdp_merge_and_extract_components</code> . The CI in question here for the total number of mutations assigned to an averaged raw cluster.
<code>ground.truth.sig</code>	Optional. Either a string with the path to file with ground truth signatures or and ICAMS catalog with the ground truth signatures. These are the signatures used to construct the ground truth spectra.
<code>ground.truth.exp</code>	Optional. Ground truth exposure matrix or path to file with ground truth exposures. If NULL skip checks that need this information.
<code>overwrite</code>	If TRUE overwrite <code>out.dir</code> if it exists, otherwise raise an error.
<code>out.dir</code>	Directory that will be created for the output; if <code>overwrite</code> is FALSE then abort if <code>out.dir</code> already exists.
<code>gamma.alpha</code>	Shape parameter of the gamma distribution prior for the Dirichlet process concentration parameters; in this function the gamma distributions for all Dirichlet processes, except possibly the top level process, are the same.
<code>gamma.beta</code>	Inverse scale parameter (rate parameter) of the gamma distribution prior for the Dirichlet process concentration parameters; in this function the gamma distributions for all Dirichlet processes, except possibly the top level process, are the same. We recommend <code>gamma.alpha = 1</code> and <code>gamma.beta = 20</code> for single-base-substitution signatures extraction; <code>gamma.alpha = 1</code> and <code>gamma.beta = 50</code> for doublet-base-substitution/INDEL signature extraction
<code>gamma0.alpha</code>	See figure B.1 from Nicola Robert's thesis. The shape parameter (α_0) of the gamma distribution priors for the Dirichlet process concentration parameters (γ_0) for G_0 .

<code>gamma0.beta</code>	See figure B.1 from Nicola Robert's thesis. Inverse scale parameter (rate parameter, β_0) of the gamma distribution priors for the Dirichlet process concentration parameters (γ_0) for G_0 .
<code>checkpoint.chlist</code>	If TRUE, checkpoint the (unclean) chlist to "initial.chlist.Rdata" in the current working directory. and checkpoint the clean chlist to "clean.chlist.Rdata" in the current working directory.
<code>checkpoint.l.chain</code>	If TRUE checkpoint the sample chain to current working directory, in a file called <code>sample.chain.seed_number.Rdata</code> .
<code>prior.sigs</code>	A matrix containing prior signatures.
<code>prior.pseudoc</code>	A numeric list. Pseudo counts of each prior signature. Recommended is 1000. In practice, it may be advisable to put lower weights on prior signatures that you do not expect to be present in your dataset, or even exclude some priors entirely.
<code>burnin.multiplier</code>	A checkpoint setting. <code>burnin.multiplier</code> rounds of burnin iterations will be run. After each round, a burn-in chain will be save for checkpoint. A total number of 10,000 iterations is recommended for most analysis. However, number of iterations can be adjusted based on the size of dataset. According to our experience, 20,000 iterations are needed when analyzing all PCAWG7 genomes (2,780 samples). The burnin can be continued from a checkpoint file with ExtendBurnin .
<code>burnin.checkpoint</code>	Default is False. If True, a checkpoint for burnin will be created.

Value

Invisibly, a list with the following elements:

- signature** The extracted signature profiles as a matrix; rows are mutation types, columns are samples (e.g. tumors).
- exposure** The inferred exposures as a matrix of mutation counts; rows are signatures, columns are samples (e.g. tumors).
- multi.chains** A `hdpSampleMulti-class` object. This object has the method `chains` which returns a list of `hdpSampleChain-class` objects. Each of these sample chains objects has a method `final_hdpState` (actually the methods seems to be just `hdp`) that returns the `hdpState` from which it was generated.
- sum_raw_clusters_after_cos_merge** A matrix containing aggregated spectra of raw clusters after cosine similarity merge step in [hdp_merge_and_extract_components](#).
- sum_raw_clusters_after_nonzero_categ** A matrix containing aggregated spectra of raw clusters after non-zero category selecting step in [hdp_merge_and_extract_components](#).
- clust_hdp0_ccc4** A matrix containing aggregated spectra of raw clusters moving to `hdp.0` after non-zero category selection step in [hdp_merge_and_extract_components](#).
- clust_hdp0_ccc5** A matrix containing aggregated spectra of raw clusters moving to `hdp.0` after non-zero observation selection step in [hdp_merge_and_extract_components](#).

PlotSamplesHighSigExp

Plot hdp signature exposure in each sample. This function returns the plot of top 5 samples with the highest exposure to a signature. Each spectrum's title is in the form of: SampleName(Proportion of Signature Assginment) This function is here because it is specific for signature extraction application.

Description

Plot hdp signature exposure in each sample. This function returns the plot of top 5 samples with the highest exposure to a signature. Each spectrum's title is in the form of: SampleName(Proportion of Signature Assginment) This function is here because it is specific for signature extraction application.

Usage

```
PlotSamplesHighSigExp(
  retval,
  hdpsample,
  input.catalog,
  col_comp = NULL,
  incl_numdata_plot = TRUE,
  ylab_numdata = "Number of data items",
  ylab_exp = "Component exposure",
  leg.title = "Component",
  cex.names = 0.6,
  cex.axis = 0.7,
  mar = c(1, 4, 2, 0.5),
  oma = c(1.5, 1.5, 1, 1)
)
```

Arguments

retval	An object return from extract_ccc_cdc_from_hdp
hdpsample	A hdpSampleChain-class or hdpSampleMulti-class object including output from extract_components_from_clusters
input.catalog	Input spectra catalog as a matrix or in ICAMS format.
col_comp	Colours of each component, from 0 to the max number. If NULL, default colors will be used
incl_numdata_plot	Logical - should an upper barplot indicating the number of data items per DP be included? (Default TRUE)
ylab_numdata	Vertical axis label for numdata plot
ylab_exp	Vertical exis label for exposure plot
leg.title	Legend title
cex.names	Expansion factor for bar labels (dpnames) in exposure plot
cex.axis	Expansion factor for vertical-axis annotation

mar See ?par
oma See ?par

PrepInit	<i>Initialize hdp object Allocate process index for hdp initialization. Prepare for hdp_init</i>
----------	--

Description

Initialize hdp object Allocate process index for hdp initialization. Prepare for [hdp_init](#)

Usage

```
PrepInit(
  multi.types,
  input.catalog,
  verbose = TRUE,
  K.guess,
  gamma.alpha = 1,
  gamma.beta = 1,
  gamma0.alpha = gamma.alpha,
  gamma0.beta = gamma.beta
)
```

Arguments

<code>multi.types</code>	<p>A logical scalar or a character vector. If FALSE, The HDP analysis will regard all input spectra as one tumor type. HDP structure as one parent node for all tumors</p> <p>If TRUE, the HDP analysis will infer tumor types based on the string before ":" in their names. e.g. tumor type for "Ovary-AdenoCA::S.500" would be "Ovary-AdenoCA" HDP structure as a grandparent node for whole data and one parent node for each tumor type</p> <p>If <code>multi.types</code> is a character vector, then it should be of the same length as the number of columns in <code>input.catalog</code>, and each value is the name of the tumor type of the corresponding column in <code>input.catalog</code>.</p> <p>If not FALSE, HDP will give a parent node for each tumor type and a grandparent node for the whole dataset.</p> <p>e.g. <code>c("Ovary-AdenoCA", "Kidney-RCC")</code>.</p>
<code>input.catalog</code>	Input spectra catalog as a matrix or in ICAMS format.
<code>verbose</code>	If TRUE then message progress information.
<code>K.guess</code>	Suggested initial value of the number of clusters. Usually, the number of clusters is two times of the number of extracted signatures. Passed to dp_activate as <code>initcc</code> .
<code>gamma.alpha</code>	Shape parameter of the gamma distribution prior for the Dirichlet process concentration parameters; in this function the gamma distributions for all Dirichlet processes, except possibly the top level process, are the same.

<code>gamma.beta</code>	<p>Inverse scale parameter (rate parameter) of the gamma distribution prior for the Dirichlet process concentration parameters; in this function the gamma distributions for all Dirichlet processes, except possibly the top level process, are the same.</p> <p>We recommend <code>gamma.alpha = 1</code> and <code>gamma.beta = 20</code> for single-base-substitution signatures extraction; <code>gamma.alpha = 1</code> and <code>gamma.beta = 50</code> for doublet-base-substitution/INDEL signature extraction</p>
<code>gamma0.alpha</code>	See figure B.1 from Nicola Robert's thesis. The shape parameter (α_0) of the gamma distribution priors for the Dirichlet process concentration parameters (γ_0) for G_0 .
<code>gamma0.beta</code>	See figure B.1 from Nicola Robert's thesis. Inverse scale parameter (rate parameter, β_0) of the gamma distribution priors for the Dirichlet process concentration parameters (γ_0) for G_0 .

PriorSetupAndActivate

Generate an HDP Gibbs sampling chain from a spectra catalog.

Description

Generate an HDP Gibbs sampling chain from a spectra catalog.

Usage

```
PriorSetupAndActivate(
  prior.sigs,
  prior.pseudoc,
  gamma.alpha = 1,
  gamma.beta = 1,
  K.guess,
  gamma0.alpha = gamma.alpha,
  gamma0.beta = gamma.beta,
  multi.types = F,
  input.catalog,
  verbose = TRUE,
  seedNumber = 1
)
```

Arguments

<code>prior.sigs</code>	A matrix containing prior signatures.
<code>prior.pseudoc</code>	A numeric list. Pseudo counts of each prior signature. Recommended is 1000. In practice, it may be advisable to put lower weights on prior signatures that you do not expect to be present in your dataset, or even exclude some priors entirely.
<code>gamma.alpha</code>	Shape parameter of the gamma distribution prior for the Dirichlet process concentration parameters; in this function the gamma distributions for all Dirichlet processes, except possibly the top level process, are the same.

<code>gamma.beta</code>	Inverse scale parameter (rate parameter) of the gamma distribution prior for the Dirichlet process concentration parameters; in this function the gamma distributions for all Dirichlet processes, except possibly the top level process, are the same.
<code>K.guess</code>	Suggested initial value of the number of signatures, passed to <code>dp_activate</code> as <code>initcc</code> .
<code>gamma0.alpha</code>	See figure B.1 from Nicola Robert's thesis. The shape parameter (α_0) of the gamma distribution priors for the Dirichlet process concentration parameters (γ_0) for G_0 .
<code>gamma0.beta</code>	See figure B.1 from Nicola Robert's thesis. Inverse scale parameter (rate parameter, β_0) of the gamma distribution priors for the Dirichlet process concentration parameters (γ_0) for G_0 .
<code>multi.types</code>	A logical scalar or a character vector. If <code>FALSE</code> , The HDP analysis will regard all input spectra as one tumor type. If <code>TRUE</code> , the HDP analysis will infer tumor types based on the string before "::" in their names. e.g. tumor type for "SA.Syn.Ovary-AdenoCA::S.500" would be "SA.Syn.Ovary-AdenoCA" If <code>multi.types</code> is a character vector, then it should be of the same length as the number of columns in <code>input.catalog</code> , and each value is the name of the tumor type of the corresponding column in <code>input.catalog</code> . e.g. <code>c("SA.Syn.Ovary-AdenoCA", "SA.Syn.Kidney-RCC")</code> .
<code>input.catalog</code>	Input spectra catalog as a matrix or in <code>ICAMS</code> format.
<code>verbose</code>	If <code>TRUE</code> then message progress information.
<code>seedNumber</code>	A random seeds passed to <code>dp_activate</code> .

Value

Invisibly, an `hdpState-class` object as returned from `dp_activate`.

RunHdpxParallel	<i>Extract mutational signatures and optionally compare them to existing signatures and exposures.</i>
-----------------	--

Description

Extract mutational signatures and optionally compare them to existing signatures and exposures.

Usage

```
RunHdpxParallel(
  input.catalog,
  seedNumber = 123,
  K.guess,
  multi.types = FALSE,
  verbose = TRUE,
  burnin = 5000,
  burnin.multiplier = 2,
  burnin.checkpoint = FALSE,
```

```

post.n = 200,
post.space = 100,
post.cputer = 3,
post.verbosity = 0,
CPU.cores = 20,
num.child.process = 20,
cos.merge = 0.9,
confident.prop = 0.9,
noise.prop = 0.5,
hc.cutoff = 0.1,
ground.truth.sig = NULL,
ground.truth.exp = NULL,
overwrite = TRUE,
out.dir = NULL,
gamma.alpha = 1,
gamma.beta = 20,
gamma0.alpha = gamma.alpha,
gamma0.beta = gamma.beta,
checkpoint.chlist = TRUE,
checkpoint.l.chain = TRUE,
prior.sigs = NULL,
prior.pseudoc = NULL
)

```

Arguments

input.catalog	Input spectra catalog as a matrix or in ICAMS format.
seedNumber	A random seeds passed to dp_activate .
K.guess	Suggested initial value of the number of clusters. Usually, the number of clusters is two times of the number of extracted signatures. Passed to dp_activate as <code>initcc</code> .
multi.types	<p>A logical scalar or a character vector. If <code>FALSE</code>, The HDP analysis will regard all input spectra as one tumor type. HDP structure as one parent node for all tumors</p> <p>If <code>TRUE</code>, the HDP analysis will infer tumor types based on the string before "::" in their names. e.g. tumor type for "Ovary-AdenoCA::S.500" would be "Ovary-AdenoCA" HDP structure as a grandparent node for whole data and one parent node for each tumor type</p> <p>If <code>multi.types</code> is a character vector, then it should be of the same length as the number of columns in <code>input.catalog</code>, and each value is the name of the tumor type of the corresponding column in <code>input.catalog</code>.</p> <p>If not <code>FALSE</code>, HDP will give a parent node for each tumor type and a grandparent node for the whole dataset.</p> <p>e.g. <code>c("Ovary-AdenoCA", "Kidney-RCC")</code>.</p>
verbose	If <code>TRUE</code> then message progress information.
burnin	Pass to hdp_burnin burnin. The number of burn-in iterations
burnin.multiplier	A checkpoint setting. <code>burnin.multiplier</code> rounds of burnin iterations will be run. After each round, a burn-in chain will be save for checkpoint. A total number of 10,000 iterations is recommended for most analysis. However,

	number of iterations can be adjusted based on the size of dataset. According to our experience, 20,000 iterations are needed when analyzing all PCAWG7 genomes (2,780 samples). The burnin can be continued from a checkpoint file with ExtendBurnin .
<code>burnin.checkpoint</code>	Default is False. If True, a checkpoint for burnin will be created.
<code>post.n</code>	Pass to hdp_posterior_sample <code>n</code> . The number of posterior samples to collect.
<code>post.space</code>	Pass to hdp_posterior_sample <code>space</code> . The number of iterations between collected samples.
<code>post.cptiter</code>	Pass to hdp_posterior_sample and hdp_burnin <code>cptiter</code> . The number of iterations of concentration parameter sampling to perform after each iteration
<code>post.verbosity</code>	Pass to hdp_posterior_sample <code>verbosity</code> . Verbosity of debugging statements. No need to change unless for development purpose
<code>CPU.cores</code>	Number of CPUs to use; this should be no more than <code>num.child.process</code> .
<code>num.child.process</code>	Number of posterior sampling chains; can set to 1 for testing. We recommend 20 for real data analysis
<code>cos.merge</code>	The cosine similarity threshold for merging raw clusters from the posterior sampling chains into "components" i.e. signatures; passed to extract_components_from_clusters
<code>confident.prop</code>	Pass to interpret_components . clusters with at least <code>confident.prop</code> of posterior samples are high confident signatures
<code>noise.prop</code>	Pass to interpret_components . Clusters with less than <code>noise.prop</code> of posterior samples are noise signatures
<code>hc.cutoff</code>	Pass to extract_components_from_clusters . The cutoff of height of hierarchical clustering dendrogram
<code>ground.truth.sig</code>	Optional. Either a string with the path to file with ground truth signatures or and ICAMS catalog with the ground truth signatures. These are the signatures used to construct the ground truth spectra.
<code>ground.truth.exp</code>	Optional. Ground truth exposure matrix or path to file with ground truth exposures. If NULL skip checks that need this information.
<code>overwrite</code>	If TRUE overwrite <code>out.dir</code> if it exists, otherwise raise an error.
<code>out.dir</code>	Directory that will be created for the output; if <code>overwrite</code> is FALSE then abort if <code>out.dir</code> already exists.
<code>gamma.alpha</code>	Shape parameter of the gamma distribution prior for the Dirichlet process concentration parameters; in this function the gamma distributions for all Dirichlet processes, except possibly the top level process, are the same.
<code>gamma.beta</code>	Inverse scale parameter (rate parameter) of the gamma distribution prior for the Dirichlet process concentration parameters; in this function the gamma distributions for all Dirichlet processes, except possibly the top level process, are the same. We recommend <code>gamma.alpha = 1</code> and <code>gamma.beta = 20</code> for single-base-substitution signatures extraction; <code>gamma.alpha = 1</code> and <code>gamma.beta = 50</code> for doublet-base-substitution/INDEL signature extraction

<code>gamma0.alpha</code>	See figure B.1 from Nicola Robert's thesis. The shape parameter (α_0) of the gamma distribution priors for the Dirichlet process concentration parameters (γ_0) for G_0 .
<code>gamma0.beta</code>	See figure B.1 from Nicola Robert's thesis. Inverse scale parameter (rate parameter, β_0) of the gamma distribution priors for the Dirichlet process concentration parameters (γ_0) for G_0 .
<code>checkpoint.chlist</code>	If TRUE, checkpoint the (unclean) chlist to "initial.chlist.Rdata" in the current working directory. and checkpoint the clean chlist to "clean.chlist.Rdata" in the current working directory.
<code>checkpoint.l.chain</code>	If TRUE checkpoint the sample chain to current working directory, in a file called <code>sample.chain.seed_number.Rdata</code> .
<code>prior.sigs</code>	A matrix containing prior signatures.
<code>prior.pseudoc</code>	A numeric list. Pesudo counts of each prior signature. Recommended is 1000. In practice, it may be advisable to put lower weights on prior signatures that you do not expect to be present in your dataset, or even exclude some priors entirely.

Value

Invisibly, a list with the following elements:

- signature** The extracted signature profiles as a matrix; rows are mutation types, columns are signatures including high confident signatures - 'hdp' signatures and moderate confident signatures - 'potential hdp' signatures.
- signature.post.samp.number** A dataframe with two columns. The first column corresponds to each signature in `signature` and the second columns contains the number of posterior samples that found the raw clusters contributing to the signature.
- signature.cdc** A `comp_dp_counts` like dataframe. Each column corresponds to the sum of all `comp_dp_counts` matrices of the raw clusters contributing to each signature in `codesignature`
- exposureProbs** The inferred exposures as a matrix of mutation probabilities; rows are signatures, columns are samples (e.g. tumors).
- noise.signature** The extracted signature profiles as a matrix; rows are mutation types, columns are signatures with less than `noise.prop` of posterior samples
- noise.post.samp.number** A data frame with two columns. The first column corresponds to each signature in `noise.signature` and the second column contains the number of posterior samples that found the raw clusters contributing to the signature.
- noise.cdc** A `comp_dp_counts` like data frame. Each column corresponds to the sum of all `comp_dp_counts` matrices of the raw clusters contributing to each signature in `code-noise.signature`
- extracted.retval** A list object returned from `codeinterpret_components`.

SetupAndActivate *Generate an HDP Gibbs sampling chain from a spectra catalog.*

Description

Generate an HDP Gibbs sampling chain from a spectra catalog.

Usage

```
SetupAndActivate(
  input.catalog,
  seedNumber = 1,
  K.guess,
  multi.types = FALSE,
  verbose = TRUE,
  gamma.alpha = 1,
  gamma.beta = 1,
  gamma0.alpha = gamma.alpha,
  gamma0.beta = gamma.beta
)
```

Arguments

<code>input.catalog</code>	Input spectra catalog as a matrix or in ICAMS format.
<code>seedNumber</code>	A random seeds passed to dp_activate .
<code>K.guess</code>	Suggested initial value of the number of clusters. Usually, the number of clusters is two times of the number of extracted signatures. Passed to dp_activate as <code>initcc</code> .
<code>multi.types</code>	<p>A logical scalar or a character vector. If <code>FALSE</code>, The HDP analysis will regard all input spectra as one tumor type. HDP structure as one parent node for all tumors</p> <p>If <code>TRUE</code>, the HDP analysis will infer tumor types based on the string before "::" in their names. e.g. tumor type for "Ovary-AdenoCA::S.500" would be "Ovary-AdenoCA" HDP structure as a grandparent node for whole data and one parent node for each tumor type</p> <p>If <code>multi.types</code> is a character vector, then it should be of the same length as the number of columns in <code>input.catalog</code>, and each value is the name of the tumor type of the corresponding column in <code>input.catalog</code>.</p> <p>If not <code>FALSE</code>, HDP will give a parent node for each tumor type and a grandparent node for the whole dataset.</p> <p>e.g. <code>c("Ovary-AdenoCA", "Kidney-RCC")</code>.</p>
<code>verbose</code>	If <code>TRUE</code> then message progress information.
<code>gamma.alpha</code>	Shape parameter of the gamma distribution prior for the Dirichlet process concentration parameters; in this function the gamma distributions for all Dirichlet processes, except possibly the top level process, are the same.

<code>gamma.beta</code>	<p>Inverse scale parameter (rate parameter) of the gamma distribution prior for the Dirichlet process concentration parameters; in this function the gamma distributions for all Dirichlet processes, except possibly the top level process, are the same.</p> <p>We recommend <code>gamma.alpha = 1</code> and <code>gamma.beta = 20</code> for single-base-substitution signatures extraction; <code>gamma.alpha = 1</code> and <code>gamma.beta = 50</code> for doublet-base-substitution/INDEL signature extraction</p>
<code>gamma0.alpha</code>	See figure B.1 from Nicola Robert's thesis. The shape parameter (α_0) of the gamma distribution priors for the Dirichlet process concentration parameters (γ_0) for G_0 .
<code>gamma0.beta</code>	See figure B.1 from Nicola Robert's thesis. Inverse scale parameter (rate parameter, β_0) of the gamma distribution priors for the Dirichlet process concentration parameters (γ_0) for G_0 .

Value

Invisibly, an `hdpState-class` object as returned from `dp_activate`.

`SetupAndPosterior` *Generate an HDP Gibbs sampling chain from a spectra catalog.*

Description

Generate an HDP Gibbs sampling chain from a spectra catalog.

Usage

```
SetupAndPosterior(
  input.catalog,
  seedNumber = 1,
  K.guess,
  multi.types = FALSE,
  verbose = TRUE,
  burnin = 5000,
  post.n = 50,
  post.space = 50,
  cpiter = 3,
  post.verbosity = 0,
  gamma.alpha = 1,
  gamma.beta = 20,
  gamma0.alpha = gamma.alpha,
  gamma0.beta = gamma.beta,
  checkpoint.1.chain = TRUE,
  burnin.multiplier = 2,
  burnin.checkpoint = TRUE,
  prior.sigs = NULL,
  prior.pseudoc = NULL
)
```

Arguments

<code>input.catalog</code>	Input spectra catalog as a matrix or in ICAMS format.
<code>seedNumber</code>	A random seeds passed to dp_activate .
<code>K.guess</code>	Suggested initial value of the number of clusters. Usually, the number of clusters is two times of the number of extracted signatures. Passed to dp_activate as <code>initcc</code> .
<code>multi.types</code>	<p>A logical scalar or a character vector. If <code>FALSE</code>, The HDP analysis will regard all input spectra as one tumor type. HDP structure as one parent node for all tumors</p> <p>If <code>TRUE</code>, the HDP analysis will infer tumor types based on the string before ":" in their names. e.g. tumor type for "Ovary-AdenoCA::S.500" would be "Ovary-AdenoCA" HDP structure as a grandparent node for whole data and one parent node for each tumor type</p> <p>If <code>multi.types</code> is a character vector, then it should be of the same length as the number of columns in <code>input.catalog</code>, and each value is the name of the tumor type of the corresponding column in <code>input.catalog</code>.</p> <p>If not <code>FALSE</code>, HDP will give a parent node for each tumor type and a grandparent node for the whole dataset.</p> <p>e.g. <code>c("Ovary-AdenoCA", "Kidney-RCC")</code>.</p>
<code>verbose</code>	If <code>TRUE</code> then message progress information.
<code>burnin</code>	Pass to hdp_burnin <code>burnin</code> . The number of burn-in iterations
<code>post.n</code>	Pass to hdp_posterior_sample <code>n</code> . The number of posterior samples to collect.
<code>post.space</code>	Pass to hdp_posterior_sample <code>space</code> . The number of iterations between collected samples.
<code>cpiter</code>	Pass to hdp_burnin <code>cpiter</code> . The number of iterations of concentration parameter sampling to perform after each iteration.
<code>post.verbosity</code>	Pass to hdp_posterior_sample <code>verbosity</code> . Verbosity of debugging statements. No need to change unless for development purpose
<code>gamma.alpha</code>	Shape parameter of the gamma distribution prior for the Dirichlet process concentration parameters; in this function the gamma distributions for all Dirichlet processes, except possibly the top level process, are the same.
<code>gamma.beta</code>	<p>Inverse scale parameter (rate parameter) of the gamma distribution prior for the Dirichlet process concentration parameters; in this function the gamma distributions for all Dirichlet processes, except possibly the top level process, are the same.</p> <p>We recommend <code>gamma.alpha = 1</code> and <code>gamma.beta = 20</code> for single-base-substitution signatures extraction; <code>gamma.alpha = 1</code> and <code>gamma.beta = 50</code> for doublet-base-substitution/INDEL signature extraction</p>
<code>gamma0.alpha</code>	See figure B.1 from Nicola Robert's thesis. The shape parameter (α_0) of the gamma distribution priors for the Dirichlet process concentration parameters (γ_0) for G_0 .
<code>gamma0.beta</code>	See figure B.1 from Nicola Robert's thesis. Inverse scale parameter (rate parameter, β_0) of the gamma distribution priors for the Dirichlet process concentration parameters (γ_0) for G_0 .

checkpoint.1.chain	If TRUE checkpoint the sample chain to current working directory, in a file called <code>sample.chain.seed_number.Rdata</code> .
burnin.multiplier	A checkpoint setting. <code>burnin.multiplier</code> rounds of burnin iterations will be run. After each round, a burn-in chain will be save for checkpoint. A total number of 10,000 iterations is recommended for most analysis. However, number of iterations can be adjusted based on the size of dataset. According to our experience, 20,000 iterations are needed when analyzing all PCAWG7 genomes (2,780 samples). The burnin can be continued from a checkpoint file with ExtendBurnin .
burnin.checkpoint	Default is False. If True, a checkpoint for burnin will be created.
prior.sigs	A matrix containing prior signatures.
prior.pseudoc	A numeric list. Pesudo counts of each prior signature. Recommended is 1000. In practice, it may be advisable to put lower weights on prior signatures that you do not expect to be present in your dataset, or even exclude some priors entirely.
post.cpiter	Pass to hdp_posterior_sample and hdp_burnin <code>cpiter</code> . The number of iterations of concentration parameter sampling to perform after each iteration

Value

Invisibly, an [hdpSampleChain-class](#) object as returned from [hdp_posterior](#).

```
test.ground.truth.sig
      test.ground.truth.sig
```

Description

Ground truth signature for testing.
Ground truth exposure for testing.

Usage

```
test.ground.truth.sig
test.ground.truth.exposure
```

Format

An [ICAMS](#) catalog (each column is a COSMICv3 signature, each row is a mutation type, e.g. ACT > AGT).

An exposure matrix (each row is a COSMICv3 signature used to generate `test.spectra`, each column is a sample).

test.spectra	<i>test.spectra</i>
--------------	---------------------

Description

Synthetic SBS96 spectra for testing.

Usage

```
test.spectra
```

Format

An [ICAMS](#) catalog (each column is a sample, each row is a mutation type, e.g. ACT > AGT).

Index

* datasets

test.ground.truth.sig, 27
test.spectra, 28

AnalyzeAndPlotretval, 2

ChainBurnin, 3
chains, 4, 8, 16
ChainsDiagnosticPlot, 4
CleanChlist, 4
CombineChainsAndExtractSigs, 5, 8
CombinePosteriorChains, 6
comp_dp_counts, 6, 23
ComponentDiagnosticPlotting, 8

dp_activate, 3, 12, 14, 18, 20, 21, 24–26

ExtendBurnin, 3, 9, 12, 16, 22, 27
extract_ccc_cdc_from_hdp, 17
extract_components_from_clusters,
6, 17, 22

final_hdpState, 4, 8, 16

GenerateAverageCluster, 9
Generateppindex, 10
GeneratePriorppindex, 10

hdp_burnin, 3, 9, 12, 15, 21, 22, 26, 27
hdp_extract_components, 7, 15
hdp_init, 18
hdp_merge_and_extract_components,
7, 8, 15, 16
hdp_posterior, 9, 27
hdp_posterior_sample, 12, 14, 15, 22,
26, 27
hdpState-class, 3, 9

ICAMS, 2, 5, 7, 8, 10–12, 14, 15, 17, 18,
20–22, 24, 26–28
interpret_components, 6, 22, 23

MultipleSetupAndPosterior, 5, 11

OldRunHdpParallel, 13

PlotSamplesHighSigExp, 17
PrepInit, 18
PriorSetupAndActivate, 19

RunHdpParallel, 13, 20

SetupAndActivate, 24
SetupAndPosterior, 25

test.ground.truth.exposure
(test.ground.truth.sig), 27
test.ground.truth.sig, 27
test.spectra, 28